

Digital Cell Morphology by CellaVision

CellaVision® Proficiency Software



The CellaVision Proficiency Software is an innovative tool designed to help laboratory managers assess, monitor and promote staff competency in the area of cell morphology. This web-based tool makes it easier to set up a structured and effective process for proficiency assessment and competency promotion without requiring too much time, effort or resources.

When implemented on its own or together with CellaVision hardware and supporting applications, the CellaVision Proficiency Software acts as a platform for competency promotion in the area of cell morphology, facilitating ambitions and efforts to deliver **more accurate** and **standardized** results.

Digital Cell Morphology by CellaVision

CellaVision® Proficiency Software

How does it work?

- The software is designed to evaluate staff members ability to perform blood and body fluid differentials
- The web-based tool can be accessed from any computer connected to the internet, using secure log-in to personalized user accounts
- The examiner can create test scenarios by uploading digital slides, either from their own CellaVision DM database or from a library of digital slides supplied by CellaVision
- The participants are automatically notified when a new test is available and instructed to perform differentials on one or more digital slides
- Upon completion, participants can browse their test result, compare their differential to that of the examiner and to those done by their co-participants
- The examiner is automatically provided with a performance report which highlights areas for future training and competency development

Why CellaVision Proficiency Software?

- Suitable for any size organization
- Can be tailored to explore identified competency gaps or meet particular training needs and goals
- Tests are quick and easy to create and manage
- Tests are quick and easy for participants to access and complete
- Results are automatically compiled
- Easy to revisit previous tests in order to keep track of participants progress

Features:

- Web-based, with instant web browser access from anywhere
- Integrated notification mechanism to participants
- Automated reporting to examiner
- Automated feedback provision to participants

FOR A MORE DETAILED PRESENTATION, PLEASE VISIT WWW.CELLAVISION.COM/PRODUCTS

TECHNICAL SPECIFICATIONS:

RECOMMENDED PC SPECIFICATIONS

- Operating systems:
Win Vista / Win 7 / Win 8 / Win 10
- Browser:
Latest version of Internet Explorer
- 2 GB RAM
- 8 Mbps bandwidth
- Screen resolution 1920x1080

DEPENDENCIES

- Microsoft Silverlight
- Transfer Tool (can only be used with CellaVision Remote Review Software Version 3.1 or higher)

EDITIONS

- Professional Edition (designed for most laboratories)
- Enterprise Edition (integrated health networks and EQA organizations)
- Trial Edition (free demo version)