

Meanwhile in the world of morphology ...



Ceaseless.

The CellaVision DM9600 is a true powerhouse for automated and digital workflow in cell morphology. The loading capacity is 96 slides with continuous feed and the throughput is approximately 30 slides per hour. The system automates the repetitive and sometimes messy tasks:

- scans and reads the slide barcode
- queries the Laboratory Information System
- applies oil on the slide
- places the slide automatically under the microscope
- starts the analysis

Apart from a little oil the DM9600 requires minimal attention. The system works ceaselessly regardless of volumes and working hours, freeing up time for the Medical Technologist and improving overall turnaround time.

Cease less in your workflow with DM9600!

[Read More](#)

CELLAVISION

www.cellavision.com

Look out for the next issue in mid-September.