



RAL Stainer MCDh™

Methanol free

As hematology laboratories become increasingly automated the limitations of traditional staining methods become more and more clear.

Discover MCDh™ staining, a unique and innovative formulation of Romanowsky stains created by CellaVision – RAL Diagnostics.

In combination with the RAL Stainer, Kit RAL Stainer MCDh™ allows a reproducible and standardized staining for your lab. RAL Stainer is a standalone stainer, that is both compact and efficient. Smart reagent management systems enable labs to strictly follow lab procedures. RAL Stainer is the ideal flexible solution for laboratories to automate their staining of peripheral blood smears.


KEY FEATURES

Ready-to-use


Fast

Safe


Digital quality




MCDh staining with RAL Stainer of peripheral blood smear from patient with normal CBC, analysis with DC-1 instrument



Eosinophil granulocyte



Lymphocyte



Neutrophil granulocyte

Kit RAL Stainer MCDh®

AVAILABLE KITS



Kit RAL Stainer MCDh™
 MCDh 1
 MCDh 2
 MCDh 3
 MCDh 4, concentrated

REF.: 360200-0000
 1 x 200 mL
 2 x 200 mL
 1 x 200 mL
 1 x 1L

TECHNICAL SPECIFICATIONS

APPLICATION AREAS

- Hematology

NUMBER OF TESTS

- 500 slides or 21 days

NUMBER OF BOTTLES

- 5 bottles

PRODUCT DESCRIPTION

- MCDh 1: Clear dark blue solution
REF.: 313590
- MCDh 2: Clear colourless solution
REF.: 3135702A
REF.: 3135703A
- MCDh 3: Clear dark blue solution
REF.: 313560
- MCDh 4, concentrated: Clear colourless solution
REF.: 313610

EQUIPMENT

- Works only with RAL Stainer